

Iluminador para lupas estereoscópicas



Luz Fría con 60Led Ultra Bríllate
LED Ring Light (60 LED Bulb, For Stereo Microscope)
Input voltage: 90-265V, 50/60Hz
Output: 24V 6W
Color temperature: 6400K
Lamp Life: Minimum 10,000 hours
Operating temperature: 32-104 F(0-40C)
Humidity range: 0-80%
Max. Opening: 62mm
Fuse:1A
The Power BOX dimension: 103x61x32mm
Brightness is adjustable

<https://www.youtube.com/watch?v=AgqOfBD2lvM>

{youtube}AgqOfBD2lvM{/youtube}
